

SEALFOAM®

Absorbable Polysaccharide Hemostat



Nature's Hemostat

SealFoam® Features

- 100% Plant-Based
- Fast Hemostasis
- Rapid Absorption in the Body
- Naturally Adhesive - Easy to Apply
- Non-pyrogenic, Latex Free

Clinical Benefits

- Eliminates the risks associated with human/animal-derived hemostats
- Effectively controls oozing and bleeding
- Minimizes the occurrence of re-bleeding
- Reduces operation time

Green Technology




SealFoam[®] Polysaccharide Hemostat

SealFoam[®] is a breakthrough hemostatic pad composed of AMP[®] material, a patented technology that modifies purified plant starch into ultra-hydrophilic, adhesive forming hemostatic polymers. AMP[®] material has a molecular structure that rapidly absorbs water from blood, causing a high concentration of platelets, red blood cells and coagulation proteins (thrombin, fibrinogen etc.). This dehydration process accelerates the normal clotting cascade and SealFoam[®] forms a gelled, adhesive matrix that provides a mechanical barrier to further bleeding. The ultra-hydrophilic and adhesive properties of SealFoam[®] are uniquely formulated to control bleeding and minimize the possibility of re-bleeding. Patient safety is enhanced since SealFoam[®] contains no human or animal components and is rapidly absorbed in the body within a few days. SealFoam[®] is CE marked and registered in major international markets.


SealFoam[®] Application Technique



1. **DRY:** Gauze dry the wound and remove all excess blood.



2. **APPLY:** Immediately after removing excess blood, apply SealFoam[®] over the entire wound using forceps or dry gloved fingers.



3. **COMPRESS:** For moderate bleeding, apply direct pressure over SealFoam[®] for 1 to 2 minutes.



4. **IRRIGATE:** SealFoam[®] after achieving hemostasis. Optional.

SealFoam[®] Ordering Information

Model	Reference No	Specifications (mm)	Packaging
SealFoam [®] HD	PD644	60x40x4	10 units/box
SealFoam [®] HD	PF642	60x40x2	10 units/box

SealFoam[®] is not approved for distribution in the U.S. SealFoam[®] is manufactured under patents pending in the U.S. and other countries.



2150 Ringwood Ave. San Jose | CA 95131 USA
tel: 408-428-9818 | fax: 408-383-9189
email: info@starchmedical.com