

## APPLICATION TECHNIQUE

1

### DRY

Gauze dry the wound and remove all excess blood.

2

### APPLY

Immediately after removing excess blood, apply **SUPERCLOT®** over the entire wound.

3

### COMPRESS

Apply direct pressure over **SUPERCLOT®** using moist gauze for 1 to 2 minutes

4

### IRRIGATE

Remove excess **SUPERCLOT®** after hemostasis is achieved.



## ORDER INFORMATION

Reference No.	Specifications	Packaging
SC0003	3g	5 pcs / box
SC0005	5g	5 pcs / box
SA01	200 mm	5 pcs / box
SA02	380 mm	5 pcs / box

### PATENTS

**US PATENT # US 9687501 B2**

**EU PATENT # EP 2 203 053 B1**

**JP PATENT # 5883895**

**CHINA PATENT # ZL 200810033239.3**

**INDIAN PATENT # 293613**

 **STARCH MEDICAL**  
Hemostatic Solutions

2150 Ringwood Ave.  
San Jose CA 95131 USA  
Tel: 408-428-9818  
Fax: 408-383-9189  
www.starchmedical.com

Distributor

CE<sub>1023</sub>

# SUPERCLOT®

ABSORBABLE POLYSACCHARIDE HEMOSTAT

From the leader in powder hemostats  
comes the latest innovation  
on the market.

## WHAT IS SUPERCLOT®?

### SUPERCLOT® ABSORBABLE POLYSACCHARIDE HEMOSTAT

is a medical device composed of absorbable modified polymer (AMP®) particles and delivery applicator. AMP® particles are biocompatible, non-pyrogenic and derived from purified plant starch. The device contains no human or animal components. **SUPERCLOT®** is intended as an absorbable hemostat system to control bleeding during surgical procedures or following traumatic injuries. Latest studies demonstrate that starch based hemostatic powder like **SUPERCLOT®** also **reduce postoperative adhesions.**



## STARCH BASED HEMOSTAT




## HOW DOES SUPERCLOT® WORK?

AMP® particles have a molecular structure that rapidly absorbs water from the blood. This dehydration process causes a high concentration of platelets, red blood cells, and coagulation proteins (thrombin, fibrinogen, etc.) which accelerates the normal, physiologic clotting cascade.

In contact with blood, AMP® particles support the formation of a gelled, adhesive matrix which provides a mechanical barrier to further control bleeding. Absorption occurs within a few days as AMP® particles are degraded by amylase and glucoamylase.

### GEL FORMING ADHESIVENESS DATA COMPARING SUPERCLOT® WITH OTHER HEMOSTATIC POWDERS

Gel forming adhesiveness was tested at 25% of each powder's maximum water absorption capacity.



\* **Conclusion:** SUPERCLOT® is 5 times more adhesive than the competitors